



Original

100% Made in Germany

LIGAMED® Heel Shoe

CE art. no. 32210-001

**Protection without perspiration,
nothing works better:**

- featherweight - airy - comfortable
- easy and quick handling
- anatomically accurate for all feet
- hygienic, rejection of dirt
- boil-proof, sterilisable
- re-usable
- maintenance-free


A real saver of time and money!






Original 100% Made in Germany LIGAMED® Heel Shoe

CE art. no. 32210-001



Its inner life consists of **LIGASANO® green**, a internationally in medicine approved dressings cover of a very strong polyamide grid fabric. The size of the heel shoe is 27.5 x 9 x 18 cm.



The original **LIGAMED® Heel Shoe** is nearly completely permeable to air and light as a feather. No accumulation of heat and moisture, you hardly feel to wear it.



The three-dimensional padding is anatomically accurate. The result is an optimal protection, no dents, a perfect fit (also in case of anomalies), high wearing comfort at almost every size and form of feet.



The fastening of the Heel Shoe is quick and easy. The fixation takes place in seconds with two extra long Velcro® fasteners and can adjust at will.

The original LIGAMED® Heel Shoe - safe all-round protection for every feet

place of application	field of application	protective functions	comfort	hygiene
<ul style="list-style-type: none"> • clinic and care • operating table • medical check-up • intensiv care • nursing homes • nursing services • general care • bed, wheel chair, etc. • patient transport ambulance • at home • long-term thermal stability from -20 °C to +60 °C 	<ul style="list-style-type: none"> • fixation of wound dressings • accompanied wound treatment • prevention • pressure ulcer at heels • ankles • ideal also for diabetics • no risks in use and contraindications known • LIGAMED® Heel Shoe is not a shoe for walking! 	<ul style="list-style-type: none"> • protection against pressure and knock for the whole foot in lying and seating positions • protection against fungal infections and moist feet • no hyperthermia • reduction of shearing forces 	<ul style="list-style-type: none"> • completely permeable to air • no accumulation of heat and moisture • featherweight • anatomic-like formed • fits big and small feet • high wearing comfort • maintenance-free 	<ul style="list-style-type: none"> • rejection of dirt, drying very fast • washable up to 95 °C with common household detergents • disinfection by hot wash • sterilisable with steam 4 min. at 134 °C • disinfect with common household disinfectants, alcohol- and chlorine-free • no medium for fungi and bacteria



Application tip: the LIGAMED® Heel Shoe and LIGASANO® white with a thickness of 1cm are an ideal combination. By this you use the therapeutical effect of LIGASANO® white and the Heel Shoe remains clean for a long time. That saves you time and money.